

BioSync — Live Demo Reading Cheat Sheet

Type these in manually during the demo • One level at a time

■ NORMAL — No Alert

Heart Rate	72 BPM	Blood Pressure Sys	118
Blood Pressure Dia	76	SpO2	97 %
Temperature	98.4 °F	Glucose	95 mg/dL
HRV	55 ms		

No notifications should fire. Dashboard stays green.

■ CAUTION — Mild Alert

Heart Rate	104 BPM	Blood Pressure Sys	128
Blood Pressure Dia	82	SpO2	94 %
Temperature	99.6 °F	Glucose	115 mg/dL
HRV	38 ms		

Caution-level notifications fire. Yellow/amber badges appear.

■ WARNING — Moderate Alert

Heart Rate	118 BPM	Blood Pressure Sys	142
Blood Pressure Dia	91	SpO2	91 %
Temperature	101.2 °F	Glucose	195 mg/dL
HRV	29 ms		

Warning notifications fire. Orange badges. Caregiver alerted.

■ DANGER — Critical Alert

Heart Rate	158 BPM	Blood Pressure Sys	183
Blood Pressure Dia	122	SpO2	87 %
Temperature	103.4 °F	Glucose	322 mg/dL
HRV	19 ms		

DANGER alerts fire across all metrics. Red badges. Emergency notifications sent.